

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody HM918

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	HM918
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	HBsAg
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	Borate Buffer, pH 9.2
Specificity	HBsAg, subtypes ad & ay
Matched Pairs	HM917
Preservatives	0.09% Sodium Azide