

Product Datasheet

Alkaline Phosphatase (ALP) Antibody HM920

Product Name	Alkaline Phosphatase (ALP) Antibody
Catalog Number	HM920
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Bovine Alkaline Phosphatase
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	Phosphate Buffered Saline, pH 7.4
Specificity	Alkaline Phosphatase
Preservatives	0.09% Sodium Azide