

Product Datasheet

Digoxigenin (DIG) Antibody HM921

Product Name	Digoxigenin (DIG) Antibody
Catalog Number	HM921
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Digoxigenin-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Digoxigenin
Preservatives	0.09% Sodium Azide