

Product Datasheet

Dengue Virus (DENV-pan NS1) Antibody HM931

Product Name	Dengue Virus (DENV-pan NS1) Antibody
Catalog Number	HM931
Species/Host	CHO Cells
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 with 150 mM NaCl
Specificity	NS1 protein of all 4 serotypes
Matched Pairs	HM513 & HM932 (capture)
Preservatives	0.1% Proclin 300