

Product Datasheet

Saint Louis Encephalitis Virus (SLEV-NS1) Antibody HM940

Product Name	Saint Louis Encephalitis Virus (SLEV-NS1) Antibody
Catalog Number	HM940
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant SLEV NS1
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Japanese Encephalitis Virus NS1 and Tick-Borne Encephalitis Virus NS1
Specificity	Saint Louis Encephalitis Virus
Matched Pairs	HM941, HM942 or HM943 (detection)
Preservatives	0.05% Sodium Azide