

Product Datasheet

Anti-h Ferritin 8803 SPRN-5 100111

Product Name	Anti-h Ferritin 8803 SPRN-5
Catalog Number	100111
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG2b
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Associated Products	Native Ferritin (Human Liver) antigen 270-40
Buffer	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
IEF Profile	6.4-7.5
Cross Reactivity	Human ferritin (liver) 167 %, Human ferritin (placenta) 201 %, Human ferritin (heart) < 1.4 %
Specificity	Antibody recognizes human ferritin
Shelf Life	24 months