

## Product Datasheet

### Adenovirus Antibody HM948

<b>Product Name</b>	Adenovirus Antibody
<b>Catalog Number</b>	HM948
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Purified Native Adenovirus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with Influenza A, Influenza B, RSV, Rotavirus, Astrovirus, Norovirus GI/GII, Enterovirus
<b>Specificity</b>	Adenovirus
<b>Matched Pairs</b>	HM948 (self-pairs)
<b>Preservatives</b>	0.09% Sodium Azide