

Product Datasheet

Clostridium difficile toxin B Antibody HM954

Product Name	Clostridium difficile toxin B Antibody
Catalog Number	HM954
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified native C. difficile toxin B
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with C. difficile toxin A
Specificity	C. difficile toxin B
Matched Pairs	HM953 (capture)
Preservatives	0.09% Sodium Azide