

Product Datasheet

Dengue Virus (DENV-pan NS1) Antibody HM960

Product Name	Dengue Virus (DENV-pan NS1) Antibody
Catalog Number	HM960
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Recombinant Dengue 3 NS1
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with recombinant NS1 Japanese Encephalitis Virus, West Nile Virus, Yellow Fever Virus. Minimal reactivity with recombinant NS1 Zika Virus.
Specificity	NS1 protein of all 4 serotypes
Matched Pairs	HM959 or HM961 (capture)
Preservatives	0.05% Sodium Azide