

Product Datasheet

Dengue Virus (DENV-pan NS1) Antibody HM962

| | |
|----------------------------|---------------------------------------------------------------------------------------------------|
| Product Name | Dengue Virus (DENV-pan NS1) Antibody |
| Catalog Number | HM962 |
| Species/Host | Mouse |
| Tested Applications | ELISA |
| Immunogen | Recombinant Dengue 3 NS1 |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | 2-8°C |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Purity | >95% |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Cross Reactivity | Minimal reactivity with recombinant NS1 Japanese Encephalitis Virus, West Nile Virus, Zika Virus. |
| Specificity | NS1 protein of all 4 serotypes |
| Matched Pairs | HM959 (capture) |
| Preservatives | 0.05% Sodium Azide |