

## Product Datasheet

### Human IL-12p70 ELISA Set 851.550

<b>Product Name</b>	Human IL-12p70 ELISA Set
<b>Catalog Number</b>	851.550
<b>Description</b>	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
<b>Range</b>	6.25 -200 pg/ml
<b>Sensitivity</b>	2.2 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	IL-12
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines
<b>Specificity</b>	Recognizes both natural and recombinant human IL-12. Specific to p70 heterodimer.
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant