

Product Datasheet

Parainfluenza Virus type 1 (PIV-1) Antibody J117

Product Name	Parainfluenza Virus type 1 (PIV-1) Antibody
Catalog Number	J117
Species/Host	Goat
Tested Applications	ELISA, WB, IHC, IFA
Immunogen	Parainfluenza Virus 1 Strain Cantell
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with Parainfluenza Virus 2 & 3
Specificity	Parainfluenza Virus 1 structural proteins
Preservatives	0.1% Sodium Azide