

Product Datasheet

Rotavirus Antibody J314

Product Name	Rotavirus Antibody
Catalog Number	J314
Species/Host	Mouse
Tested Applications	LF, ELISA, Immunostaining
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with Influenza A, Influenza B, RSV, Parainfluenza 1, 2 & 3, Adenovirus, M. pneumoniae
Specificity	Strains RRV, WA & bovine
Matched Pairs	J340 (detection)
Preservatives	0.1% Sodium Azide