

Product Datasheet

Cryptosporidium parvum Antibody J350

Product Name Cryptosporidium parvum Antibody

Catalog Number J350

Species/Host Mouse

Tested Applications LF, ELISA, IFA

Brand EastCoast Bio

Form/Appearance Protein A Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG1

Clonality Monoclonal

Purity >90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Specificity Intact C. parvum oocysts

Preservatives 0.1% Sodium Azide