

Product Datasheet

Escherichia coli (E. coli) O157 Antibody J352

Product Name	Escherichia coli (E. coli) O157 Antibody
Catalog Number	J352
Species/Host	Mouse
Tested Applications	LF, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG3
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	E.coli O157 (intact or extracted)
Matched Pairs	HM016 (for lateral flow assay)
Preservatives	0.1% Sodium Azide