

## Product Datasheet

### Escherichia coli (E. coli) O157 Antibody J352

<b>Product Name</b>	Escherichia coli (E. coli) O157 Antibody
<b>Catalog Number</b>	J352
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, IFA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG3
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Specificity</b>	E.coli O157 (intact or extracted)
<b>Matched Pairs</b>	HM016 (for lateral flow assay)
<b>Preservatives</b>	0.1% Sodium Azide