

Product Datasheet

Neisseria meningitidis Antibody JP025

Product Name	Neisseria meningitidis Antibody
Catalog Number	JP025
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Neisseria meningitidis whole cells (types A, B, C)
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, purified IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	May react with related microbes (antiserum in unabsorbed)
Specificity	All Neisseria Meningitidis Antigens (Types A, B, C)
Preservatives	0.1% Sodium Azide