

## Product Datasheet

### Streptococcus pneumoniae Antibody JP034

<b>Product Name</b>	Streptococcus pneumoniae Antibody
<b>Catalog Number</b>	JP034
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	LF, ELISA, WB, IHC, IFA
<b>Immunogen</b>	Whole cell mix of multiple serotypes
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Cross Reactivity</b>	May react with related microbes (antiserum in unabsorbed)
<b>Specificity</b>	Known Reactivity With S. Pneumoniae Types 3,4,6,7,9,14,18,19 & 23
<b>Preservatives</b>	0.1% Sodium Azide