

Product Datasheet

Mouse Immunoglobulin (IgG) Antibody JP038

Product Name	Mouse Immunoglobulin (IgG) Antibody
Catalog Number	JP038
Species/Host	Goat
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Immunoaffinity purified Polyclonal Antibody
Storage	2-8°C
Clonality	Polyclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Preservatives	0.02% Sodium Azide