

---

## Product Datasheet

### Rabbit Immunoglobulin (IgG) Antibody JP039

<b>Product Name</b>	Rabbit Immunoglobulin (IgG) Antibody
<b>Catalog Number</b>	JP039
<b>Species/Host</b>	Goat
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Immunoaffinity Purified Polyclonal Antibody
<b>Storage</b>	2-8°C
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline , pH 7.4
<b>Preservatives</b>	0.02% Sodium Azide