

## Product Datasheet

### Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody-Biotin JP040

<b>Product Name</b>	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody-Biotin
<b>Catalog Number</b>	JP040
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Native HIV antigen
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody-Biotin Conjugate
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4 with 50% glycerol
<b>Matched Pairs</b>	HM083 or HM084