

Product Datasheet

Bovine Serum Albumin (BSA) Antibody JP041

Product Name	Bovine Serum Albumin (BSA) Antibody
Catalog Number	JP041
Species/Host	Rabbit
Tested Applications	ELISA, IP, WB
Immunogen	Highly Purified Bovine Serum Albumin
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, delipidated serum
Storage	2-8°C. Do not freeze.
Clonality	Polyclonal
Specificity	Single precipitin arc using immunoelectrophoresis against BSA and the albumin in bovine serum
Preservatives	0.09% Sodium Azide