

## Product Datasheet

### Bovine Serum Albumin (BSA) Antibody JP042

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Bovine Serum Albumin (BSA) Antibody   |
| <b>Catalog Number</b>      | JP042   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Tested Applications</b> | ELISA, WB, IHC  |
| <b>Immunogen</b>           | Highly Purified Bovine Serum Albumin  |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Purified Polyclonal Antibody  |
| <b>Storage</b>             | 2-8°C. Do not freeze.   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Specificity</b>         | Single precipitin arc using immunoelectrophoresis against BSA and the albumin in bovine serum |
| <b>Preservatives</b>       | 0.09% Sodium Azide  |