

Product Datasheet

Ractopamine (RAC) Antibody JP045

Product Name	Ractopamine (RAC) Antibody
Catalog Number	JP045
Species/Host	Rabbit
Tested Applications	LF, ELISA
Immunogen	Ractopamine-BSA
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	EIA titer 1:500 - 1:5000.
Clonality	Polyclonal
Preservatives	0.1% Sodium Azide