

Product Datasheet

Zolpidem Antibody JP046

Product Name Zolpidem Antibody

Catalog Number JP046

Species/Host Rabbit

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Protein A Purified Polyclonal Antibody, Ig Fraction

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Clonality Polyclonal

Buffer 0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride

Specificity Zolpidem

Matched Pairs LA155 (capture)

Preservatives 0.05% Sodium Azide