

## Product Datasheet

### Zolpidem Antibody JP048

**Product Name** Zolpidem Antibody

**Catalog Number** JP048

**Species/Host** Rabbit

**Tested Applications** LF, ELISA

**Brand** EastCoast Bio

**Form/Appearance** Protein A Purified Polyclonal Antibody, Ig Fraction

**Storage** Short Term ( $\leq 2$  weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

**Clonality** Polyclonal

**Buffer** 0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride

**Specificity** Zolpidem

**Matched Pairs** LA155 (capture)