

Product Datasheet

Heartworm, Canine Antibody JP051

Product Name	Heartworm, Canine Antibody
Catalog Number	JP051
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	R. G. antigen
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	2-8°C
Clonality	Polyclonal
Buffer	0.05M Phosphate Buffered Saline, pH 7.4
Matched Pairs	FA306
Preservatives	0.1% Sodium Azide