

## Product Datasheet

### C-Reactive Protein, Canine (cCRP) Antibody JP057

<b>Product Name</b>	C-Reactive Protein, Canine (cCRP) Antibody
<b>Catalog Number</b>	JP057
<b>Species/Host</b>	Chicken
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Canine CRP isolated from dog serum
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified complete IgY fraction
<b>Storage</b>	2-8°C
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.1M Sodium Phosphate Buffer, pH 7.2 with 150 mM Sodium Chloride
<b>Preservatives</b>	0.09% Sodium Azide