

---

## Product Datasheet

### Chicken Antibody JP059

<b>Product Name</b>	Chicken Antibody
<b>Catalog Number</b>	JP059
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, WB, ICC
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Polyclonal Antibody
<b>Storage</b>	2-8°C. Do not freeze.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Preservatives</b>	0.02% Sodium Azide