

Product Datasheet

Progesterone-17-OH Antibody JP060

Product Name	Progesterone-17-OH Antibody
Catalog Number	JP060
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	17-OH Progesterone-3-CMO-BSA
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody; neat serum
Storage	≤-20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Cross Reactivity	17-OH Progesterone-3 100%, 17 α OH Progesterone < 0.1%, Progesterone < 0.1%
Preservatives	None