

## Product Datasheet

### Human Growth Hormone (hGH) Antibody JP062

<b>Product Name</b>	Human Growth Hormone (hGH) Antibody
<b>Catalog Number</b>	JP062
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Human Growth Hormone
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Polyclonal Antibody; neat serum
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Cross Reactivity</b>	HGH 100%, hPrI $\leq 0.3\%$ , hFSH $\leq 0.2\%$ , hLH $\leq 0.2\%$ , hTSH $\leq 0.2\%$
<b>Specificity</b>	Human
<b>Preservatives</b>	None