

Product Datasheet

Cortisol-3 Antibody JP064

Product Name	Cortisol-3 Antibody
Catalog Number	JP064
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Cortisol-3-BSA
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody; neat serum
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Cross Reactivity	Cortisol 100%, Prednisolone $\leq 36.0\%$, 11-Desoxycortisol $\leq 5.7\%$, Corticosterone $\leq 3.3\%$
Preservatives	None