

Product Datasheet

Cortisol-3 Antibody JP064

| | |
|----------------------------|--|
| Product Name | Cortisol-3 Antibody |
| Catalog Number | JP064 |
| Species/Host | Rabbit |
| Tested Applications | ELISA |
| Immunogen | Cortisol-3-BSA |
| Brand | EastCoast Bio |
| Form/Appearance | Polyclonal Antibody; neat serum |
| Storage | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Clonality | Polyclonal |
| Cross Reactivity | Cortisol 100%, Prednisolone $\leq 36.0\%$, 11-Desoxycortisol $\leq 5.7\%$, Corticosterone $\leq 3.3\%$ |
| Preservatives | None |