

Product Datasheet

Aldosterone-3 Antibody JP065

Product Name	Aldosterone-3 Antibody
Catalog Number	JP065
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Aldosterone-3-BSA
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, Neat Serum
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Cross Reactivity	11-Deoxycorticosterone $\leq 1.1\%$
Preservatives	None