

Product Datasheet

Streptococcus Group A Antibody JP086

Product Name	Streptococcus Group A Antibody
Catalog Number	JP086
Species/Host	Rabbit
Tested Applications	LF, ELISA
Immunogen	S. pyogenes
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	2-8°C
Clonality	Polyclonal
Purity	≥95%
Buffer	0.02M Sodium Phosphate, 0.15M Sodium Chloride
Cross Reactivity	No reactivity with tested strains of Streptococcus Group B, C, F, G, S. aureus, H. influenzae, N. perflava
Specificity	Group A Streptococcus
Preservatives	0.05% Sodium Azide