

## Product Datasheet

### Alpha 1 Acid Glycoprotein (A1GP) Antibody JP087

<b>Product Name</b>	Alpha 1 Acid Glycoprotein (A1GP) Antibody
<b>Catalog Number</b>	JP087
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, IEP, TIA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Polyclonal Antibody, IgG Fraction
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	$\geq 95\%$
<b>Buffer</b>	20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl
<b>Specificity</b>	Human A1GP
<b>Preservatives</b>	0.05% Sodium Azide