

Product Datasheet

Alpha 1 Acid Glycoprotein (A1GP) Antibody JP087

Product Name	Alpha 1 Acid Glycoprotein (A1GP) Antibody
Catalog Number	JP087
Species/Host	Goat
Tested Applications	ELISA, IEP, TIA
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, IgG Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	$\geq 95\%$
Buffer	20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl
Specificity	Human A1GP
Preservatives	0.05% Sodium Azide