

Product Datasheet

Ceruloplasmin Antibody JP088

Product Name Ceruloplasmin Antibody

Catalog Number JP088

Species/Host Goat

Tested Applications ELISA, IEP, TIA

Brand EastCoast Bio

Form/Appearance Polyclonal Antibody, IgG Fraction

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Clonality Polyclonal

Buffer 20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl

Specificity Human Ceruloplasmin

Preservatives 0.05% Sodium Azide