

Product Datasheet

Cystatin C Antibody JP089

Product Name Cystatin C Antibody

Catalog Number JP089

Species/Host Rabbit

Tested Applications ELISA

Brand EastCoast Bio

Form/Appearance Liquid; Purified Polyclonal Antibody

Storage Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Clonality Polyclonal

Purity >90%

Buffer 10 mM Phosphate Buffered Saline, pH 7.4