

Product Datasheet

6-Monoacetylmorphine (6-MAM) Antibody JP093

Product Name	6-Monoacetylmorphine (6-MAM) Antibody
Catalog Number	JP093
Species/Host	Rabbit
Tested Applications	LF, ELISA
Immunogen	6-Monoacetylmorphine-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Cross Reactivity	Diacetylmorphine <2.5%
Specificity	6-Monoacetylmorphine
Matched Pairs	170031 (capture)
Preservatives	0.05% Sodium Azide