

## Product Datasheet

### Human Immunoglobulin E (IgE) Antibody JP094

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Human Immunoglobulin E (IgE) Antibody  |
| <b>Catalog Number</b>      | JP094  |
| <b>Species/Host</b>        | Goat   |
| <b>Tested Applications</b> | ELISA, WB, IHC   |
| <b>Immunogen</b>           | Purified Human IgE Heavy ( $\epsilon$ ) chain                                      |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Affinity Purified Polyclonal Antibody  |
| <b>Storage</b>             | 2-8°C  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Purity</b>              | $\geq 95\%$ by SDS-PAGE  |
| <b>Buffer</b>              | 10 mM Sodium Phosphate, pH 7.2 with 0.15 M Sodium Chloride                         |
| <b>Cross Reactivity</b>    | No reactivity with Human immunoglobulin light chains or other Human serum proteins |
| <b>Specificity</b>         | Human IgE Heavy ( $\epsilon$ ) chain   |
| <b>Preservatives</b>       | 0.05% Sodium Azide   |