

## Product Datasheet

### Human Kappa Light Chain Antibody-AP JP095

<b>Product Name</b>	Human Kappa Light Chain Antibody-AP
<b>Catalog Number</b>	JP095
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, WB, IHC, IF
<b>Immunogen</b>	Purified Human Kappa ( $\kappa$ ) Light chain
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Polyclonal Antibody
<b>Storage</b>	2-8°C
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥95% by SDS-PAGE
<b>Buffer</b>	10 mM Sodium Phosphate, pH 7.2 with 0.15 M Sodium Chloride
<b>Cross Reactivity</b>	No reactivity with Human immunoglobulin heavy chains, Human immunoglobulin lambda ( $\lambda$ ) light chains or other Human serum proteins
<b>Specificity</b>	Human Kappa ( $\kappa$ ) Light chains
<b>Preservatives</b>	0.05% Sodium Azide