

Product Datasheet

Human Immunoglobulin D (IgD) Antibody JP097

Product Name	Human Immunoglobulin D (IgD) Antibody
Catalog Number	JP097
Species/Host	Goat
Tested Applications	ELISA, FC
Immunogen	Human IgD
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	2-8°C
Isotype	IgG
Clonality	Polyclonal
Buffer	Borate Buffered Saline, pH 8.2
Cross Reactivity	Cross Adsorbed with Human IgA, IgM, IgG. Possible reactivity with IgD from other species.
Specificity	Human IgD
Preservatives	None