

## Product Datasheet

### Human Immunoglobulin D (IgD) Antibody JP097

<b>Product Name</b>	Human Immunoglobulin D (IgD) Antibody
<b>Catalog Number</b>	JP097
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, FC
<b>Immunogen</b>	Human IgD
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Polyclonal Antibody
<b>Storage</b>	2-8°C
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	Borate Buffered Saline, pH 8.2
<b>Cross Reactivity</b>	Cross Adsorbed with Human IgA, IgM, IgG. Possible reactivity with IgD from other species.
<b>Specificity</b>	Human IgD
<b>Preservatives</b>	None