

Product Datasheet

K2 Antibody JP098

| | |
|----------------------------|--|
| Product Name | K2 Antibody |
| Catalog Number | JP098 |
| Species/Host | Rabbit |
| Tested Applications | LF, ELISA |
| Immunogen | K2-BgG |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Polyclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Clonality | Polyclonal |
| Buffer | 0.05 M Phosphate Buffer, pH 7.15 with 0.85% Sodium Chloride |
| Cross Reactivity | JWH-018 100%; AM-2232 35.7%; JWH-073 10.1%; JWH-200 3.7%; JWH-019 3.3; JWH-203 <0.2%, JWH-147, JWH-175, HU-210, CP-47,947, AH-7921 <0.1% |
| Specificity | K2 |
| Matched Pairs | LA288 (capture) |
| Preservatives | 0.05% Sodium Azide |