

## Product Datasheet

### K2 Antibody JP099

<b>Product Name</b>	K2 Antibody
<b>Catalog Number</b>	JP099
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	K2-BgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Polyclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.05 M Phosphate Buffer, pH 7.15 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	JWH-018 100%; AM-2232 167%; JWH-073 35.7%; JWH-200 24.1%; JWH-019 5.0%; JWH-203, JWH-147, JWH-175, HU-210, CP-47,947, AH-7921 <0.1%
<b>Specificity</b>	K2
<b>Matched Pairs</b>	LA288 (capture)
<b>Preservatives</b>	0.05% Sodium Azide