

Product Datasheet

K2 Antibody JP100

Product Name	K2 Antibody
Catalog Number	JP100
Species/Host	Rabbit
Tested Applications	LF, ELISA
Immunogen	K2-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05 M Phosphate Buffer, pH 7.15 with 0.85% Sodium Chloride
Cross Reactivity	JWH-018 100%; AM-2232 12.3%; JWH-200 12.0%; JWH-073 9.1%; JWH-019 6.1; JWH-203 <0.2%; JWH-147, JWH-175, HU-210, CP-47,947, AH-7921 <0.1%
Specificity	K2
Matched Pairs	LA288 (capture)
Preservatives	0.05% Sodium Azide