

## Product Datasheet

### Hydrocodone Antibody JP102

<b>Product Name</b>	Hydrocodone Antibody
<b>Catalog Number</b>	JP102
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Hydrocodone-BgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Polyclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.05 M Phosphate Buffer, pH 7.15 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	Hydrocodone 100%; Hydromorphone 11.8%; Codeine <0.4%; Dihydrocodeine <0.2%; Morphine, Oxycodone <0.1%
<b>Specificity</b>	Hydrocodone
<b>Matched Pairs</b>	LA289 (capture)
<b>Preservatives</b>	0.05% Sodium Azide