
Product Datasheet

Albumin Antibody JP106

Product Name	Albumin Antibody
Catalog Number	JP106
Species/Host	Goat
Tested Applications	ELISA, IEP, TIA
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	2-8°C
Clonality	Polyclonal
Buffer	20 mM Sodium Phosphate, 150 mM NaCl, pH 7.2
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.05% Sodium Azide