

---

## Product Datasheet

### Albumin Antibody JP106

<b>Product Name</b>	Albumin Antibody
<b>Catalog Number</b>	JP106
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	ELISA, IEP, TIA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Polyclonal, Delipidated and 0.2 µm filtered
<b>Storage</b>	2-8°C
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	20 mM Sodium Phosphate, 150 mM NaCl, pH 7.2
<b>Specificity</b>	Single precipitin arc visible with immunoelectrophoresis
<b>Preservatives</b>	0.05% Sodium Azide